


Comfortex Healthcare Mattresses & Cushions

Product Listing and Indications for Use

Acute Care Comfort, Non-Mechanical Pressure Redistribution, Moderate to High Risk Prevention and Wound Care Stages I, II and Uncomplicated III's and IV's :

Rest-Q Unique High Specification design for Advanced Pressure Redistribution and Wound Care for patients weighing up to 350 lbs.

- Clinically proven to be more comfortable than Memory Foam mattresses.
- Multi-layered, Multi-Zoned for overall Pressure redistribution and support.
- Pressure Relief zones Off-Load pressure points for enhanced protection and wound treatment.
- Extra wide Firm Edging for safer bedside sitting, entering and exiting the bed.
- Advanced performance without the noise, movements or mechanical issues associated with powered systems.
- Comfortex' low-shear nylon top cover is both water and vapor proof for enhanced durability, protection from incontinence and infection control.



Rest-Q GS – GSL or GSC “Gentle Strength”

- GS designs utilize a unique blend of materials that deliver the **Gentle** protection for the smallest of patients and the **Strength** to comfortably and safely support the heaviest of patients.
- Adding special layers of HR foams provide enhanced support, balanced immersion and pressure redistribution for patients weighing less than 100 lbs. and also to those weighing more than 500 lbs.
- **Each GS mattress is made to fit and function atop any hospital bed manufacturer's frame...Better than the original.**



Rest-Q GSV High Specification design for Advanced Pressure Redistribution and Wound Care. **Uniquely shaped mattress to fit uniquely shaped bed frames...Perfectly.**

- Clinically proven to be more comfortable than Memory Foam mattresses.
- Multi-layered, Multi-Zoned for overall Pressure redistribution and support.
- Pressure Relief zones Off-Load pressure points for enhanced protection and wound treatment.
- Advanced clinical performance without the noise, movements or mechanical issues associated with powered systems.
- Extra perimeter cushioning buffers the bed-frame's edges for more comfortable EOB sitting & transferring.
- Comfortex' low-shear, non-stretch nylon top cover goes beyond fluid-proof to prevent vapor, viral or bacterial contamination.
- Heel Suspension Zone gently absorbs sensitive heels and moves with the bed's telescoping function without flexing upward.



Specialized Healthcare Surfaces

Advanced pressure redistribution and support for Bariatric care. High Resiliency Support and a non-accessible cover for safe Behavioral Health Comfort and Care.

Bariatric Care mattresses:

Rest-Q GS Plus Same design features as the standard Rest-Q and GS...plus:

- Thicker layers of High Resiliency foam expand the support capability for the largest of patients, weighing as much as 800 lbs. (8" thickness for truly Gentle Strength)
- Multi-layered, Multi-Zoned for overall Pressure redistribution and support.
- Pressure Relief zones Off-Load pressure points for enhanced protection and wound treatment.



Extended Capacity expandable mattress system

With Rest-Q GS Plus mattress

- Expandable length and width to fit expandable frames.
- Special blend of HR foam layers safely support patients weighing up to 1000 lbs.
- Tension relief Zones allow the mattress to move with the bed frame when adjusting positions.
- Quick release and store-tight extension sleeves.
- Expansion cushions can increase the sleep surface from 36" to 48". *Optional: length extension cushion.*



Behavioral Health/Memory Care mattresses:

Closed System * recognized by the Behavioral Health Facilities Design Guide

- High Resiliency foams maximize Comfort and Support.
- Available with either XtendaCare or Rest-Q core designs.
- No access construction for patient & staff safety.
- Double Reinforced Stitching maximizes strength, durability and security.



Comfortex Healthcare Mattresses & Cushions

Product Listing and Indications for Use

Basic Comfort & Pressure Redistribution:

Low Risk Prevention, Comfort and Support for persons weighing up to 350 lbs.

NewLife - Enhanced Comfort and Moderate Risk Prevention.

- 3 layers, Multi Density High Resiliency foam core.
- Softer sections beneath the shoulders and heels, while more supportive beneath the hips.
- Extra wide Firm Edging for safer bedside sitting, entering and exiting the bed.
- Comfortex' low-shear nylon top cover is water and even vapor proof for enhanced durability and incontinence protection.



Stretcher mattresses Safe, Comfortable Support

Response stretcher mattresses 4" dual density construction with an unequaled moisture proof top cover (5" thick available for **Bariatric Support**)

- Specially engineered, shear reducing nylon top cover goes beyond water proof to prevent microscopic fluid vapors from entering the mattress.
- Dual Layers of High Resiliency, High Specification foams.
- Top layer enhances comfort and pressure redistribution.
- Supportive base layer ensures extended-stay comfort and helps prevent bottoming out.
- Unique, embedded CPR boards are available.
- Built to fit. (Easy to complete order form is provided, customer's order history is retained for easy re-ordering.)




For more information, please call us at 800-445-4007

Patient Safety

Extra Wide Edging to support easier, safer, more stable patient transfers.

- 4-inch width is twice that of most other mattresses.
- The soft top layer spans over the Edging for more comfortable edge-of-bed sitting.
- Edging is bonded along the vertical & horizontal planes for far greater strength and durability.



Raised Secure Edging (SE) Ultra firm raised edging, helps redirect patients away from bed's edge and an accidental fall.

- Central Exit Zone facilitates easier, supported transferring.
- Discrete appearance preserves patient dignity.



Stretch fabric... **Why we don't.** Commonly used stretch fabrics report to be breathable, and still waterproof... ...Consider the infection control risks.

- When subjected to the weight, heat and movement of a resting patient, microscopic pores in the stretch material elongate and any vaporized fluids can be drawn into and contaminate the mattress core.
- Comfortex' non-shearing Stratos nylon DOES NOT stretch and has an unmatched 25 year history of protection & performance.
- A proven barrier against moisture, vapor and even viral pass-thru.



the **Landing Strip** bedside fall mat The Nation's first and still the Nation's best.

- High Density compressed foam core absorbs shock and prevents pass thru impact or energy rebound.
- Precision bonding and welded perimeter.
- Partially detached bottom cover maximizes energy dissipation.
- Unified, not foldable, for greater protection and durability.



For more information, please call us at 800-445-4007

Custom Responses: Comfortex is proud of its ability to design custom surfaces for unique needs. Hospitals, Health Systems, Long-Term Care Providers and the VA have all reached out to Comfortex with unique patient conditions and/or needs.

Our staff will work quickly to design & produce surfaces that can help the provider meet the unique challenge. This service is often donated to the provider in a spirit of co-operation and with an understanding that responding to the unique, is a Reflection of Excellence...An effort that improves us all.

***Collaboration can Foster Innovation...
...Innovation can Change the World.***

Just some of the unique mattresses & cushions we've developed to support our customers...support their customers

- Specialized Cancer Treatments and Care.
- Protective sleep environment for patients & families dealing with the challenges associated with Huntington's Disease.
- Protective Side-rail Seizure pads.
- Pediatric Mattresses & Cushions.



For more information, please call us at 800-445-4007